

Easy to use with minimal training

For better perception ICare

iCare IC100 tonometer

The iCare's patented rebound technology is a reliable technique for obtaining accurate and rapid IOP measurements from all your patients. Measuring requires no anesthetic drops, air puffs or other preparation.

No anesthesia, no air puff

The easy-to-use iCare IC100™ tonometer incorporates guiding features and an intuitive user interface to maximize efficiency in clinical routines. No specialized skills are needed. Just load, align, and measure. The accurate iCare IC100 tonometer is the professional's choice.

Key features and benefits



iCare AMS
Dual-mode
operation (single
or rapid 6-series
measurement)



iCare EasyNav Advanced navigation interface



OLED display
A bright OLED
display screen
guides operation
and presents
results



iCare EasyPos
Alignment indicator:
green confirms
tonometer
positioning



Functionality with comfort



Effortless probe loading using the original iCare probe.



The momentary contact of the probe with the cornea is harmless and barely noticed by the patient.



A probe cannot accidentally drop from the device.



With iCare IC100 IOP measuring is fast, accurate and reliably repeatable.



Simple user menu, available in multiple languages.



All the measurements required for calculation of a result can be performed individually or in an automated series, at the press of a button.

Technical data

iCare IC100	
Model name	TA011
Dimensions (W x H x D)	29 mm x 95 mm x 215 mm / 1.1" x 3.7" x 8.5"
Weight	135 g (without batteries) / 0.3 lbs (without batteries) 230 g (with 4 x AA batteries) / 0.5 lbs (with 4 x AA batteries)
Power supply	4 x AA non-rechargeable batteries, 1.5 V alkaline LR6
Measurement range	7 mmHg – 50 mmHg
Accuracy	± 1.2 mmHg (≤ 20 mmHg) and ± 2.2 mmHg (> 20 mmHg)
Repeatability	< 8 % (coefficient of variation)

iCare. For better perception.

iCare is a trusted partner in ophthalmic diagnostics, offering physicians fast, easy-to-use, and reliable tools for diagnosis of glaucoma, diabetic retinopathy, and macular degeneration (AMD). Our product assortment includes automated TrueColor imaging devices, perimeters and handheld rebound tonometers.

We believe that ophthalmic care should be accessible, effortless, and reliable, and we aim to establish the next level of eye care.

C€ 0598



Icare Finland Oy

Äyritie 22 01510 Vantaa, Finland Ph. +358 9 8775 1150 info@icare-world.com

Centervue S.p.A.

Via San Marco 9H 35129 Padova, Italy Ph. +39 049 501 8399 info@icare-world.com

Icare USA, Inc.

8601 Six Forks Rd., Ste 310 Raleigh, NC 27615 Ph. +1 888.422.7313 Fax +1 877.477.5485 infoUSA@icare-world.com

www.icare-world.com



iCare is a trademark of Icare Finland Oy. Icare Finland Oy, Icare USA, Inc., iCare World Australia Pty Ltd and Centervue S.p.A., are parts of Revenio Group and represent the brand iCare. iCare IC100 is a trademark or a registered trademark of Icare Finland Oy in the EU and the USA and/or other countries. Not all products, accessories, services or offers are approved in every market. Approved labelling and instructions may vary from one country to another. This material is intended for healthcare professionals only.

© 2024 Icare Finland Oy. All rights reserved.